## MAGIC NEEDLE MICRO CANNULA GAUGE CHART



## Luer Lock base for firm and greater resistance:

The base of the micro cannula needle is specially designed for a reinforced connection with the micro cannula to facilitate the use of the doctor.

The standrad luer lock base can be fitted perfectly with the cannula needle syringe.

For stable locking and enduring higher pressure.

Cannula to needle sizing				
Cannula needle size	Outer die meter	Cannula color	Matched sharp needle size	
(Blunt needle)	Outer diameter	Cannula color	(Sharp needle)	Needle color
14G X 70MM	2.0MM		20G X 25MM	
16G X 70MM	1.6MM		20G X 25MM	
18G X 70MM	1.2MM		20G X 25MM	
18G X 50MM	1.2MM		20G X 25MM	
20G X 70MM	0.9MM		22G X 25MM	
20G X 50MM	0.9MM		22G X 25MM	
21G X 70MM	0.8MM		20G X 25MM	
21G X 50MM	0.8MM		20G X 25MM	
22G X 70MM	0.7MM		20G X 25MM	
22G X 50MM	0.7MM		20G X 25MM	
23G X 70MM	0.6MM		22G X 25MM	
23G X 50MM	0.6MM		22G X 25MM	
23G X 38MM	0.6MM		22G X 25MM	
25G X 50MM	0.5MM		23G X 25MM	
25G X 38MM	0.5MM		23G X 25MM	
27G X 50MM	0.4MM		25G X 25MM	
27G X 38MM	0.4MM		25G X 25MM	
30G X 25MM	0.3MM		25G X 25MM	

Horizon (Changzhou) Medical Technology Co., Ltd. kathy@med-horizon.com www.med-horizon.com



